

SMART HOSPITAL SYSTEM

Martin Emmanuel Chang (212100199)
Moses Anthony Kwik (212100291)

Proyek - IEE3032 - Pemrograman Berbasis Web untuk IoT



TABLE OF CONTENT

**01 Overview System
dan Subsystem**

**02 Halaman Utama
dan Fiturnya**

**03 Halaman Info dan
Log1 - Log3**

**04 Halaman Admin
dan Database**

05 Machine Learning

**06 Node-Red dan
Kode MQTT**

07 Kode Django

08 Kode-kode Proyek



OVERVIEW SYSTEM

Di industri kesehatan, mengelola dan menganalisis data pasien sangat penting. Sebuah proyek Django dapat dirancang untuk mengumpulkan, menyimpan, dan menganalisis data dari berbagai sumber dalam fasilitas kesehatan. Proyek ini dapat difokuskan pada meningkatkan perawatan pasien, mengoptimalkan kondisi ruangan operasi, dan memastikan kepatuhan terhadap regulasi kesehatan.



SUBSISTEM #1 - PATIENT MONITORING SYSTEM

Location: Patient Room

Sensors: IV Level sensor (0 - 1000 mL), Heartbeat sensor (0 - 120 bpm), Fluorescence-based glucose level sensor (100 - 190 mg/dL)

Actuator: Alert, IV flow regulator, Insulin pump

Di dalam patient room, ada sistem Patient Monitoring System yang menggunakan 3 sensor dan 3 aktuator.

Cara kerja aktuator:

- Aktuator Alert akan aktif jika IV level sensor lebih rendah dari 10 mL atau Heartbeat sensor lebih rendah dari 70 bpm atau Glucose sensor lebih dari 160 mg/dL, atau semua kombinasi kondisi-kondisi ini.
- IV Flow Regulator akan aktif untuk membuat meningkatkan flow rate cairan IV jika heartbeat sensor dibawah 90 bpm dan level sensor di atas 150 mL.
- Insulin Pump akan aktif untuk meningkat insulin pump jika glucose level lebih dari 150 mg/dL.



SUBSISTEM #2 - VISITOR HEALTH CHECK SYSTEM

Location: Hospital

Sensors: Camera detection (Boolean), Microbial sensor (50 - 250 j/m³), Thermal sensor (35 - 40 C)

Actuator: Laser pointer, Door lock, Disinfectant spray

Ada sistem Visitor Health Check System yang menggunakan 3 sensor dan 3 aktuator.

- Laser pointer akan aktif jika Camera mendeteksi visitor yang tidak pakai masker
- Door lock akan aktif jika visitor tidak memakai masker, Microbial sensor mendeteksi bakteri sebanyak 150 J/m³ di visitor atau Thermal sensor mendeteksi suhu badan yang diatas 38°C, atau semua kombinasi kondisi-kondisi ini.
- Disinfectant spray akan aktif jika Microbial sensor mendeteksi bakteri sebanyak 150 J/m³ di visitor.



SUBSISTEM #3 - SURGERY ROOM CONDITIONING SYSTEM

Location: Surgery Room

Sensors: Room pressure sensor (0 - 15 Pa), Temperature sensor (20 - 30 C), Humidity sensor (0 - 80 %)

Actuator: Electric damper, HVAC system, Alarm

Ada sistem Operating Room Conditioning System yang menggunakan 3 sensor dan 3 aktuator.


- HVAC Sistem akan aktif menyalahkan cooling system jika sensor temperatur diatas 25 °C atau sensor humiditas diatas 65%, atau dua-duanya.
- Alarm akan aktif jika sensor temperatur diatas 27 °C atau sensor humiditas diatas 70% atau sensor room pressure di bawah 2.0 Pa atau semua kombinasi kondisi-kondisi ini.
- Electric damper akan aktif untuk meningkatkan tekanan ruangan jika sensor room pressure dibawah 2.5 Pa.

HALAMAN UTAMA

Smart Hospital System

INFO

Subsystem 1: Patient Monitoring System



alert_actuator IV_flow_regulator_actuator insulin_pump_actuator



ML Output: 1




ML Output: 1



ML Output: 1

[LOG](#)

Subsystem 2: Visitor Health Check System



laser_pointer_actuator door_lock_actuator disinfectant_spray_actuator



ML Output: 1




ML Output: 1



ML Output: 1

[LOG](#)

Subsystem 3: Surgery Room Conditioning System



electric_damper_actuator HVAC_system_actuator alarm_actuator



ML Output: 1



ML Output: 1



ML Output: 1

[LOG](#)

HALAMAN INFO

System Explanation

[BACK](#)

In the healthcare industry, managing and analyzing patient data is crucial. A Django project can be designed to collect, store, and analyze data from various sources within a healthcare facility. This project can be focused on improving patient care, optimizing hospital operations, and ensuring compliance with health regulations.

Patient Monitoring System

Location: Patient's room


SENSORS



IV Level Sensor: Measures the infusion fluid level in the patient's infuse bag. Threshold is 0 - 1000mL.



Heartbeat Sensor: Measures the patient's heartbeat rate. Threshold is 0 - 120 bpm.




Flourescent Glucose Based Sensor: Measures the concentration of glucose in the patient's blood stream using fluorophore. Threshold is 100 - 190 mg/dL.

Visitor Health Check System

Location: All around the hospital

SENSORS



Camera Detection Sensor: Uses facial and object recognition system to detect whether a visitor is wearing a medical face mask or not around the hospital. Threshold is 0 (mask detected) or 1 (no mask detected).



Microbial Sensor: Using UV-C light to detect presence of germs or microorganisms on a visitor before they enter the hospital. Threshold is 50 - 250 J/m³.




Thermal Sensor: Measures the visitor's body temperature before entering the hospital. Threshold is 35 - 40 °C.


Surgery Room Conditioning System

Location: Surgery room

SENSORS



Room Pressure Sensor: Measures the surgery room's air pressure. Maintaining good air pressure in a surgery room is crucial for infection control, airborne pathogen control and ventilation. Threshold is 0 - 15 Pa.






Temperature Sensor: Measures the surgery room's temperature. Maintaining good room temperature in a surgery room is crucial for patient comfort, staff comfort and equipment performance. Threshold is 20 - 30 °C.



Humidity Sensor: Measures the surgery room's humidity. Maintaining good room

HALAMAN INFO

100 - 190 mg/dL.




ACTUATORS




Alert Actuator: The alert actuator will activate if the IV level sensor value is below 10 mL or the heartbeat sensor value is below 70 bpm or the fluorescent based glucose sensor value is above 160 mg/dL, or any combination of these conditions.

IV Flow Actuator: The IV flow actuator will activate to increase the IV fluid flow rate if the heartbeat sensor value is below 90 bpm and the level sensor value is above 150 mL.

Insulin Pump Actuator: The insulin pump actuator will activate to increase insulin delivery if the glucose level is above 150 mg/dL.



temperature before entering the hospital. Threshold is 35 - 40 °C.




ACTUATORS




Laser Pointer Actuator: Laser pointer actuator is activated when the camera detection sensor returns value 1 or there is a visitor that doesn't uses a medical face mask.

Door Lock Actuator: The door lock actuator will activate if the visitor is not wearing a medical face mask, the microbial sensor detects bacteria of 150 J/m^3 on the person, or the thermal sensor detects a body temperature above 38°C, or any combination of these conditions.

Disinfectant Spray Actuator: The disinfectant spray actuator will activate if the microbial sensor detects bacteria of 150 J/m^3 on the visitor.



Humidity Sensor: Measures the room's humidity. Maintaining good room humidity in a surgery room is crucial for static electricity control, equipment performance, patient and staff comfort, and wound tissue healing. Threshold is 0 - 80 %.



ACTUATORS


Electric Damper Actuator: The electric damper actuator will activate to increase room pressure if the room pressure sensor reads below 2.5 Pa.

HVAC System Actuator: The HVAC system actuator will activate to turn on the cooling system if the temperature sensor reads above 25 °C or if the humidity sensor reads above 65%, or both.

Alarm Actuator: The alarm actuator will activate if the temperature sensor reads above 27 °C, or if the humidity sensor reads above 70%, or if the room pressure sensor reads below 2.0 Pa, or any combination of these conditions.

[BACK](#)

HALAMAN INFO



Machine Learning Technique: Logistic Regression with CSV Files

Machine Learning Technique: Logistic Regression

Logistic Regression is a statistical model used for binary classification problems. It models the probability that an instance belongs to a particular class using a logistic function.

Training Using CSV Files

CSV (Comma-Separated Values) files are used to store tabular data, with each row representing a data instance. Training involves using a labeled dataset to teach the model, and the CSV files store the features and labels for training instances.

Integration of Logistic Regression and CSV Files

Logistic regression is applied to the data stored in CSV files during the training process. The model adjusts its parameters to minimize the difference between predictions and actual labels, learning patterns from the training data.

Outcome

After training, the logistic regression model can predict the probability of an instance belonging to a class. It can be applied to new data, calculating probabilities and assigning classes based on a chosen threshold.

HALAMAN LOG1

Patient Monitoring System Log

BACK

| Timestamp | IV Level Sensor | Heartbeat Sensor | Fluorescent Based Sensor | Alert | IV Flow Regulator | Insulin Pump |
|---------------------|-----------------|------------------|--------------------------|-------|-------------------|--------------|
| 2023-12-09 22:24:14 | 143.39 | 58.0 | 152.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:24:07 | 368.27 | 74.0 | 166.0 | 0.0 | 1.0 | 0.0 |
| 2023-12-09 22:23:58 | 735.03 | 87.0 | 184.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:51 | 558.43 | 35.0 | 136.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:42 | 191.88 | 68.0 | 137.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:33 | 319.58 | 31.0 | 189.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:23 | 119.24 | 104.0 | 143.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:14 | 199.91 | 31.0 | 134.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:06 | 703.44 | 56.0 | 142.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:59 | 657.09 | 57.0 | 114.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:53 | 643.69 | 113.0 | 167.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:46 | 509.84 | 84.0 | 144.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:40 | 174.39 | 113.0 | 119.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:34 | 892.57 | 49.0 | 186.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:23 | 527.94 | 74.0 | 144.0 | 1.0 | 0.0 | 1.0 |
| 2023-12-09 22:22:21 | 63.14 | 1.0 | 125.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:22:20 | 25.46 | 110.0 | 153.0 | 0.0 | 0.0 | 1.0 |
| 2023-12-09 22:22:18 | 12.78 | 89.0 | 120.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:22:16 | 155.57 | 48.0 | 151.0 | 0.0 | 1.0 | 0.0 |
| 2023-12-09 22:22:14 | 667.51 | 12.0 | 149.0 | 0.0 | | 1.0 |
| | | | | 0.0 | | |

HALAMAN LOG2

Health Regulation System Log

BACK

| Timestamp | Camera_sensor | Microbial_sensor | Thermal_sensor | Laser_pointer | Door_lock | Disinfectant_spray |
|---------------------|---------------|------------------|----------------|---------------|-----------|--------------------|
| 2023-12-09 22:25:36 | 1.00 | 148.0 | 38.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:26 | 0.00 | 233.0 | 36.0 | 0.0 | 0.0 | 1.0 |
| 2023-12-09 22:25:15 | 0.00 | 130.0 | 35.0 | 0.0 | 0.0 | 1.0 |
| 2023-12-09 22:24:56 | 1.00 | 222.0 | 38.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:50 | 0.00 | 151.0 | 35.0 | 1.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:42 | 0.00 | 125.0 | 39.0 | 0.0 | 1.0 | 0.0 |
| 2023-12-09 22:24:39 | 0.00 | 60.0 | 37.0 | 1.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:31 | 1.00 | 233.0 | 37.0 | 1.0 | 0.0 | 1.0 |
| 2023-12-09 22:24:22 | 0.00 | 55.0 | 38.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:14 | 0.00 | 159.0 | 37.0 | 1.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:07 | 1.00 | 243.0 | 39.0 | 1.0 | 0.0 | 0.0 |
| 2023-12-09 22:23:59 | 1.00 | 170.0 | 38.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:52 | 1.00 | 229.0 | 35.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:43 | 1.00 | 220.0 | 39.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:33 | 1.00 | 169.0 | 36.0 | 0.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:23 | 0.00 | 157.0 | 35.0 | 1.0 | 0.0 | 1.0 |
| 2023-12-09 22:23:15 | 0.00 | 75.0 | 35.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:23:06 | 1.00 | 186.0 | 39.0 | 0.0 | 0.0 | 1.0 |
| 2023-12-09 22:23:00 | 0.00 | 62.0 | 35.0 | 1.0 | 0.0 | 1.0 |
| 2023-12-09 22:22:54 | 0.00 | 66.0 | 38.0 | 0.0 | 1.0 | 1.0 |


HALAMAN LOG3

Surgery Room Conditioning System Log

[BACK](#)

| Timestamp | Room_pressure_sensor | Temperature_sensor | Humidity_sensor | Electric_damper | HVAC_system | Alarm |
|---------------------|----------------------|--------------------|-----------------|-----------------|-------------|-------|
| 2023-12-09 22:26:31 | 1.33 | 23.26 | 46.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:26:24 | 13.14 | 23.19 | 50.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:26:20 | 0.94 | 22.55 | 77.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:26:18 | 6.99 | 25.12 | 43.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:26:16 | 13.51 | 24.91 | 61.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:26:14 | 11.38 | 20.50 | 38.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:26:12 | 3.14 | 27.04 | 24.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:26:09 | 9.08 | 27.00 | 39.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:26:06 | 10.16 | 22.35 | 29.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:25:58 | 13.01 | 23.50 | 11.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:50 | 10.81 | 20.04 | 0.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:43 | 14.99 | 26.74 | 31.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:36 | 3.40 | 27.07 | 21.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:27 | 0.05 | 22.02 | 5.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:25:16 | 11.09 | 20.31 | 14.0 | 1.0 | 1.0 | 1.0 |
| 2023-12-09 22:24:57 | 1.49 | 22.22 | 28.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:50 | 14.21 | 28.50 | 9.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:42 | 0.20 | 20.14 | 6.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:39 | 14.56 | 25.50 | 29.0 | 0.0 | 0.0 | 0.0 |
| 2023-12-09 22:24:32 | 1.43 | 28.92 | 23.0 | 1.0 | 1.0 | 1.0 |

HALAMAN ADMIN DAN DATABASE



HALAMAN ADMIN DAN DATABASE

Work

WebloT

Presen

Node

Site a

Selec

Chat

Mikse 2

IEE30

Lapoi

Lapoi

M15

127.0.0.1:8000/admin/SmartHospitalSystem/actuator_list/

Django administration

WELCOME, MARTIN. [VIEW SITE](#) / [CHANGE PASSWORD](#) / [LOG OUT](#)

Home > Smarthospitalsystem > Actuator_lists

Start typing to filter...

AUTHENTICATION AND AUTHORIZATION

Groups

+ Add

Users

+ Add

SMARTHOSPITALSYSTEM

Actuator_lists

+ Add

Historical actuator_lists

+ Add

Historical sensor_list_ints

+ Add

Sensor_list_ints

+ Add

Select actuator_list to change

ADD ACTUATOR_LIST +

Search

Action: -----

Go

0 of 9 selected

| <input type="checkbox"/> | NAME | DATA | TIMESTAMP |
|--------------------------|--------------------|------|---------------------|
| <input type="checkbox"/> | Alarm | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | HVAC_system | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Electric_damper | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Disinfectant_spray | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Door_lock | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Laser_pointer | 0.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | IV_flow_regulator | 1.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Insulin_pump | 1.0 | 2023-12-09 19:43:07 |
| <input type="checkbox"/> | Alert | 0.0 | 2023-12-09 19:43:07 |

9 actuator_lists

FILTER

By name

All

Alarm

Alert

Disinfectant_spray

Door_lock

Electric_damper

HVAC_system

IV_flow_regulator

Insulin_pump


Laser_pointer

By timestamp

All

2023-12-09 19:43:07

HALAMAN ADMIN DAN DATABASE



HALAMAN ADMIN DAN DATABASE

Work

WebloT

Present

Node-R

Select s

ChatGP

Mikse 2

IEE3021

Laporan

Laporan

M15 - T

127.0.0.1:8000/admin/SmartHospitalSystem/sensor_list_int/

Django administration

WELCOME, MARTIN. [VIEW SITE](#) / [CHANGE PASSWORD](#) / [LOG OUT](#)

Home > Smarthospitalsystem > Sensor_list_ints

Start typing to filter...

AUTHENTICATION AND AUTHORIZATION

Groups + Add

Users + Add

SMARTHOSPITALSYSTEM

Actuator_lists + Add

Historical actuator_lists + Add

Historical sensor_list_ints + Add

Sensor_list_ints + Add

Select sensor_list_int to change

ADD SENSOR_LIST_INT +

Search

Action: ----- Go 0 of 9 selected

| NAME | DATA | TIMESTAMP |
|--------------------------|--------------------|---------------------|
| humidity_sensor | 13.0 | 2023-12-09 19:43:07 |
| temperature_sensor | 21.347751721546267 | 2023-12-09 19:43:06 |
| room_pressure_sensor | 5.529001355587223 | 2023-12-09 19:43:06 |
| thermal_sensor | 36.0 | 2023-12-09 19:43:06 |
| Microbial_sensor | 199.0 | 2023-12-09 19:43:06 |
| camera_sensor | 0.0 | 2023-12-09 19:43:06 |
| flourescent_based_sensor | 155.0 | 2023-12-09 19:43:06 |
| heart_beat_sensor | 49.0 | 2023-12-09 19:43:06 |
| IV_level_sensor | 366.86576043584427 | 2023-12-09 19:43:06 |

9 sensor_list_ints

FILTER

By name

All

IV_level_sensor

Microbial_sensor

camera_sensor

flourescent_based_sensor

heart_beat_sensor

humidity_sensor

room_pressure_sensor

temperature_sensor

thermal_sensor


By timestamp

All


2023-12-09 19:43:06

2023-12-09 19:43:07

MACHINE LEARNING DATA



```
Project_UAS > csv > Insulin_pump.csv
1 IV_level_sensor,heart_beat_sensor,flourescent_based_sensor,Insulin_pump
2 10,89,100,0
3 49,70,120,0
4 10,95,149,0
5 10,95,150,0
6 10,95,151,0
7 20,89,200,1
8 40,80,160,1
9 50,120,198,1
10
```



MACHINE LEARNING CODE

```
for actuator in actuators_table1:
    df = pd.read_csv("S:/Uni_stuff/Semester_5/UAS WEB/Project_UAS/csv/" + str(actuator) + ".csv")

    sensor_names_table1 = [sensor.name for sensor in sensors_table1]
    features = df[sensor_names_table1]
    target = df[str(actuator)]

    X_train, X_test, y_train, y_test = train_test_split(features, target, test_size=0.2, random_s

    # Create a logistic regression model
    model = LogisticRegression()

    # Train the model
    model.fit(X_train, y_train)

    # Make predictions on the test set
    predictions = model.predict(X_test)

    # Calculate accuracy (you might want to store or log this value for evaluation)
    accuracy = accuracy_score(y_test, predictions)

    # Use the latest sensor data for prediction
    #new_data = [[latest_sensor_data_dict[sensor] for sensor in features.columns]]

    new_data = pd.DataFrame(
        np.array([[latest_sensor_data_dict[sensor] for sensor in features.columns]]).reshape(1, -1),
        columns=['IV_level_sensor', 'heart_beat_sensor', 'flourescent_based_sensor']
    )

    result = model.predict(new_data)
```


Penjelasan:

Dalam kode, nama sensor diambil dari list terus dibuat sebagai features yaitu row atas di csv data. target dipilih sebagai aktuatornya.

Model yang digunakan adalah Logistic Regression dan akan di train menggunakan data dari csv file. Kemudian model ini akan digunakan untuk predict nilai akuator dengan values sensor saat ini.



NODE-RED



MQTT KODE

```
def on_connect(client, userdata, flags, rc, ):
    print("Connected with result code " + str(rc))

    # Subscribe to topics for all sensors
    for sensor_name in sensor_names:
        client.subscribe(sensor_name)

def on_message(client, userdata, msg):
    global sensor_data_history_list

    sensor_name = msg.topic
    latest_sensor_data = float(msg.payload.decode())

    # Add the latest sensor data to the history list
    current_time = datetime.now().strftime("%Y-%m-%d %H:%M:%S")
    sensor_data_history_list.append({'sensor_name': sensor_name, 'data': latest_sensor_data, 'timestamp': current_time})

    # Keep only the last 100 sensor readings in the list
    sensor_data_history_list = sensor_data_history_list[-20:]

    # Update the latest sensor data in the dictionary
    latest_sensor_data_dict[sensor_name] = latest_sensor_data
    historical_sensor_data_dict[sensor_name].append({'data': latest_sensor_data, 'timestamp': current_time})
    historical_sensor_data_dict[sensor_name] = historical_sensor_data_dict[sensor_name][-20:]

# Set up the MQTT client
mqtt_client = mqtt.Client()
mqtt_client.on_connect = on_connect
mqtt_client.on_message = on_message
mqtt_broker_address = "127.0.0.1"
mqtt_client.connect(mqtt_broker_address, 1889)
mqtt_client.loop_start()
```



THANK YOU

